


Foodservice Equipment • Material Handling, Healthcare & Cleanroom Equipment • Retail Display

## New **QuadAdjust™** Wire Shelving



patented  
**QuadTruss®**  
design

**PATENT  
PENDING!**



channel of shelf corner fits into slot of collar attached to split sleeve



QuadAdjust™ shelf unit

Item #: \_\_\_\_\_  
Model #: \_\_\_\_\_  
Project #: \_\_\_\_\_  
SIS #: \_\_\_\_\_

Optional Accessories	Model #	Qty.
Casters		
Ledges		
Dividers		
Rods & tabs		
Aluminum split sleeves		
Shelf markers		

- Patented QuadTruss® shelves (patent #5,390,803) for increased strength and durability.
- Shelves feature special interlocking corners that fit in slots of collars mounted to split sleeves, allowing addition or removal of a shelf without the need to disassemble the entire shelf unit.
- Tapered split sleeves create a positive lock that becomes stronger as the load increases.
- Posts are numbered in increments of 1" (25mm) to ensure fast and level assembly.
- Leveling feet are provided to help compensate for uneven floor surfaces.
- Shelf utilizes a pincer-type design with the mat wire sandwiched between the two top truss wires, adding significant strength and distributing the entire load without stress and strain on the welds.


### Assembly is a simple 3-step exercise:



1) Snap the split sleeves onto a post over the number of your choice.






2) Attach the QuadAdjust™ slotted collar onto the split sleeves.



3) Slide shelf into slotted collar to complete.

Spec sheets available for viewing, printing or downloading from our online literature library at [www.eaglegrp.com](http://www.eaglegrp.com)

# QuadAdjust™ Wire Shelving



## Shelves

width x length		weight		EAGLEbrite®		chrome		Valu-Master®		Valu-Gard®		EAGLEgard®	
in.	mm	lb.	kg	model #	list	model #	list	model #	list	model #	list	model #	list
14" x 36"	356 x 914	8	3.6	QA1436Z	60	QA1436C	75	QA1436V	60	QA1436VG	60	QA1436E	84
14" x 48"	356 x 1219	10	4.6	QA1448Z	67	QA1448C	84	QA1448V	67	QA1448VG	67	QA1448E	94
14" x 60"	356 x 1524	14	6.4	QA1460Z	76	QA1460C	92	QA1460V	76	QA1460VG	76	QA1460E	103
14" x 72"	356 x 1829	17	7.7	QA1472Z	84	QA1472C	103	QA1472V	84	QA1472VG	84	QA1472E	116
18" x 36"	457 x 914	9	4.1	QA1836Z	63	QA1836C	79	QA1836V	63	QA1836VG	63	QA1836E	92
18" x 48"	457 x 1219	12	5.5	QA1848Z	71	QA1848C	89	QA1848V	71	QA1848VG	71	QA1848E	105
18" x 60"	457 x 1524	17	7.7	QA1860Z	86	QA1860C	101	QA1860V	86	QA1860VG	86	QA1860E	121
18" x 72"	457 x 1829	20	9.1	QA1872Z	95	QA1872C	115	QA1872V	95	QA1872VG	95	QA1872E	135
21" x 36"	533 x 914	11	5.0	QA2136Z	69	QA2136C	88	QA2136V	69	QA2136VG	69	QA2136E	100
21" x 48"	533 x 1219	14	6.4	QA2148Z	77	QA2148C	98	QA2148V	77	QA2148VG	77	QA2148E	114
21" x 60"	533 x 1524	18	8.2	QA2160Z	91	QA2160C	111	QA2160V	91	QA2160VG	91	QA2160E	129
21" x 72"	533 x 1829	24	10.9	QA2172Z	101	QA2172C	130	QA2172V	101	QA2172VG	101	QA2172E	153
24" x 36"	610 x 914	13	5.9	QA2436Z	72	QA2436C	92	QA2436V	72	QA2436VG	72	QA2436E	106
24" x 48"	610 x 1219	16	7.3	QA2448Z	82	QA2448C	106	QA2448V	82	QA2448VG	82	QA2448E	125
24" x 60"	610 x 1524	21	9.5	QA2460Z	94	QA2460C	120	QA2460V	94	QA2460VG	94	QA2460E	144
24" x 72"	610 x 1829	26	11.8	QA2472Z	110	QA2472C	139	QA2472V	110	QA2472VG	110	QA2472E	164

Note: MICROGARD® antimicrobial protection standard on all EAGLEgard® shelving.

## Posts

Numerically grooved in 1" (25mm) increments. Includes post cap and leveling bolt.

height	weight	EAGLEbrite®		chrome		Valu-Master®		Valu-Gard®		EAGLEgard®	
in.	mm	model	list	model	list	model	list	model	list	model	list
7"	178	P7-Z	9	P7-C	10	P7-V	11	P7-VG	11	P7-E	13
14"	356	P14-Z	11	P14-C	12	P14-V	13	P14-VG	13	P14-E	14
18"	457	P18-Z	13	P18-C	14	P18-V	15	P18-VG	15	P18-E	19
33"	838	P33-Z	14	P33-C	15	P33-V	16	P33-VG	16	P33-E	20
54"	1372	P54-Z	19	P54-C	19	P54-V	20	P54-VG	20	P54-E	21
63"	1600	P63-Z	19	P63-C	20	P63-V	21	P63-VG	21	P63-E	24
74"	1880	P74-Z	20	P74-C	22	P74-V	23	P74-VG	23	P74-E	26
86"	2184	P86-Z	25	P86-C	27	P86-V	27	P86-VG	27	P86-E	30

Notes:

- EAGLEbrite® posts are clear epoxy coated for use in dry or wet environments.
- MICROGARD® antimicrobial protection standard on all EAGLEgard® posts.
- For mobile application, add prefix "C" to model number. Example: CP14-E. See Spec Sheet EG01.05 for information about casters available.